


Certofix[®] Quinto

Maximizing venous access
for the polytrauma patient



More certainty, more safety, more success

B | BRAUN
SHARING EXPERTISE

Certofix® Quinto

Product Specifications

Central Venous Catheters

Managing polytrauma

With B. Braun's patented kink-proof guide wire, Certofix® Quinto provides 5 separate lumens to handle the complex demands of managing polytrauma. Administer drugs, give volume therapy, monitor pressure levels, conduct routine infusions, take samples – Certofix® Quinto allows the critical maximization of a central venous access site. Your key tool to support the patient in need.




Optimized section

- 1 x 12 G lumen for an effective volume therapy
- 1 x 16 G and 3 x 18 G lumens for highly differentiated intensive care therapy

Product Type	Puncture needle	Catheter diameter in F	Catheter length in cm	Lumen in G	Code No. (REF)
Quinto SB 1215	SB	12	15	D16/M,18/M ₂ /18/M ₃ /P12	4166833
Quinto V 1215	V	12	15	D16/M,18/M ₂ /18/M ₃ /P12	4166841
Quinto SB 1220	SB	12	20	D16/M,18/M ₂ /18/M ₃ /P12	4166850
Quinto V 1220	V	12	20	D16/M,18/M ₂ /18/M ₃ /P12	4166868
Quinto SB 1230	SB	12	30	D16/M,18/M ₂ /18/M ₃ /P12	4166876
Quinto V 1230	V	12	30	D16/M,18/M ₂ /18/M ₃ /P12	4166884

legend: V: Valve needle
SB: Seldinger needle and Braunula

All V-Sets available with kink-proof guide wire.



ECG-lead position monitoring

- Real-time placement verification
- On-the-spot certainty without the delays and expense of x-ray exposure
- Standard in all Certofix® catheters

- 1 A quick look at the ECG monitor shows: The P wave is clearly elevated indicating that the catheter tip is directly before the right atrium.
- 2 The catheter tip has been withdrawn 2 to 3 cm in the superior vena cava and the P wave has normalized. Now the catheter has reached its final verified position.



The new certainty

Central Venous Catheters

100% kink-proof guide wire

- 1 Flexible and firm, supple and strong, elastic and resilient ...
- 2 Tie it, twist it, curl it ...
- 3 The wire remains: 100% kink-proof, 100% intact



User benefits

- Smooth catheter insertion even in extremely difficult anatomical situations by:
 - kink-proof behavior
 - optimized tensile strength
 - ideal ratio of flexibility and resilience
- Real-time placement verification with ECG-lead technology – a standard feature

B | BRAUN
SHARING EXPERTISE

B. Braun Melsungen AG

P.O. Box 1120

D-34209 Melsungen

Tel +49(0) 56 61 71-0

www.bbraun.com

www.cvc-partner.com