

Drainobag®

Accessories

High-vacuum wound drainage systems

	Code- No.	Sales unit Pieces
Redon connection tubes 123 cm		
• double packed, sterile with variable drain connector Ch. 6–18	U217 0101	40
with variable drain connector Ch. 6–18 and slide clamp for secretion stop	U218 0101	40
with variable drain connector Ch. 6–18, slide clamp for secretion stop and lock fitting	U217 0555	40
with two variable drain connectors Ch. 6–18 on Y-piece, slide clamp for secretion stop and lock fitting	552 4938	40
with two variable drain connectors Ch. 6–18 on Y-piece	552 4946	40
Variable drain connector	552 4989	200
• Ch. 6–18 • can be used for Drainobag®/Drainofix®		
Redon-Drains		
• PVC (medical grade) • radiopaque • cross-perforated • graduated • 50 cm • double packed, sterile		
Ch. 6, 15 cm	U211 0600	100
Ch. 8,	U211 0800	
Ch. 10	U211 1000	
Ch. 12	U211 1200	
Ch. 14	U211 1400	
Ch. 16	U211 1600	
Ch. 18	U211 1800	
Redon-Drain with adaptor for variable drain connector		
• silicone • cross-perforated • 50 cm • double packed, sterile		
Ch. 18, 14 cm	552 5187	100



Drainobag® wound drainage systems
improve postoperative wound healing, eliminate the risk of
contamination and meet economical demands.

B | BRAUN
MEDICAL

B. Braun Melsungen AG
P.O. Box 1120
D-34209 Melsungen
Tel +49 (5661) 71-0
Fax +49 (5661) 71-2044


B. 03. 08. 99/1 Nr. 604 6800 Printed on chlorine-free bleached cellulose using 50% recycled paper

B | BRAUN

Drainobag®

always the right choice

High-vacuum wound drainage systems alternatively available with luer-lock or luer fitting



Luer-lock fitting for connection of secretion bottle and patient tubing

- reduces the risk of contamination
- prevents disconnection in case of tensile strain



Variable drain connector with graduation, alternatively for CH. 6-8, 10-12, 14-16 and 18.

- Safe connection of drains of any Charrière-size



Fine graduation up to 50 ml for smallest volumes of secretion in combination with a standard graduation of 10/50 ml.

- for exact determination of drained secretion volume

	Code- No.	Sales unit Pieces
Drainobag® 600 V ● with connection tube, slide clamp and variable drain connector ● double packed, sterile PVC, sterilized by ethylene oxide	0552 3605	25
Drainobag® 600 ● without connection tube PVC, sterilized by ethylene oxide	0552 3656	30
Drainobag® Lock 600 V ● with connection tube, lock connection, slide clamp and variable drain connector Ch. 6-18 ● double packed, sterile PVC-free, sterilized by irradiation PVC, sterilized by ethylene oxide	0552 2617 0552 3648	25 25
Drainobag® Lock 600 ● without connection tube ● with lock fitting PVC-free, sterilized by irradiation PVC, sterilized by ethylene oxide	0552 2668 0552 3680	30 30
Drainobag® Lock 300 V ● with connection tube, lock adaptor, slide clamp and variable drain connector, Ch. 6-18 ● double packed, sterile PVC-free, sterilized by irradiation PVC, sterilized by ethylene oxide	0552 2307 0552 2340	30
Drainobag® Lock 300 ● lock fitting ● without connection tube PVC-free, sterilized by irradiation PVC, sterilized by ethylene oxide	0552 2358 0552 2390	30
Drainobag® 300 V ● with connection tube, slide clamp and variable drain connector, Ch. 6-18 ● double packed, sterile PVC, sterilized by ethylene oxide	0552 2322	30
Drainobag® 300 ● without connection tube PVC, sterilized by ethylene oxide	0552 2374	30

	Code- No.	Sales unit Pieces
Drainobag® 150 VL ● with long connection tube, slide clamp and variable drain connector, Ch. 6-18 ● double packed, sterile PVC, sterilized by ethylene oxide	0552 3710	30
Drainobag® 150 V PVC, sterilized by ethylene oxide	0552 3702	30
Drainobag® 150 ● without connection tube PVC, sterilized by ethylene oxide	0552 3753	30
Drainobag® Lock 150 VL ● with long connection tube, lock adaptor, slide clamp and variable drain connector, Ch. 6-18 ● double packed, sterile PVC-free, sterilized by irradiation PVC, sterilized by ethylene oxide	0552 2722 0552 3737	30 30
Drainobag® Lock 150 ● with lock fitting ● without connection tube PVC-free, sterilized by irradiation PVC, sterilized by ethylene oxide	0552 2765 0552 3761	30 30
Drainobag® 150 K6 ● with connection tube, slide clamp and variable drain connector, Ch. 6-18 ● Redon drain Ch. 6 made of PVC ● single-use Redon needle Ch. 6 ● double packed, sterile PVC, sterilized by ethylene oxide	0552 3800	30
Drainobag® 150 K8 ● with connection tube, slide clamp and variable drain connector, Ch. 6-18 ● Redon drain Ch. 8 made of PVC ● single-use Redon needle Ch. 8 ● double packed, sterile PVC, sterilized by ethylene oxide	0552 3850	30