


Ureofix[®] 111

Closed urine drainage system
with bottom outlet



- Remarkable reduction of urinary tract infections due to proper function of drip-chamber, antireflux-valve and slide-clamp.
- Agreeable handling for patients and personnel due to the all-purpose bedframe hooks and the hygienic use of the bottom outlet.

Ureofix® 111

Infusion & Collection Systems

Transparent drip chamber with non-return valve

Air filters with high hydrophobic properties

- ensure the necessary pressure compensation in the drip chamber and in the urine bag

Integrated non-return valve

- prevents reflux of urine into the drip chamber and the tubing

Pasteur's way comprising drip chamber, non-return valve and slide clamp for tubing

- designed to prevent retrograde tracking of infection with catheterized patients



Suspension

Fixation bracket with high elasticity, rectangular locking position, cord and a bedside fixation plate facilitate

- single-handed use
- safe suspension to almost all types of bed frames
- multiple suspension possibilities



Bottom outlet

Its large bore internal diameter and the length of the outlet allow for:

- a quick and efficient emptying process
- free flowing of coagulants
- contamination-free single-handed operation

Presentation

Description	Code-No.	Sales unit/pcs
Ureofix® 111 closed urine drainage system with bottom outlet		
Connecting tubing/bag-volume		
100 cm 2 l	441 3229	25
120 cm 2 l	441 3237	25
100 cm 4 l	441 3326	20
100 cm 4 l with suction ball	441 3423	20
Fixation clamp for tubing Ureofix®-urine drainage systems, non-sterile	441 1404	50
Fixation plate for Ureofix® 111 and secretion bag suitable for all bedframes flexible fastening straps non-sterile	441 0602	25



Sampling port

Secured with a smooth highly flexible material which guarantees:

- multiple aspirations of urine samples
- hygienic and efficient disinfection of the surface

B | BRAUN
SHARING EXPERTISE

B. Braun Melsungen AG

P.O. Box 1120
D-34209 Melsungen
Tel +49 (56 61) 71-0
Fax +49 (56 61) 71-2777

www.bbraun.com